



**FAST • PRECISE • AFFORDABLE**  
**3D ADD-ON FOR MICROSCOPES**


- 3D Surface Analysis
- Form & Roughness Measurement
- Automatic Depth Composition
- Z Depth Measurement



*Compatible With Compound Microscopes & Stereo-Microscopes*  
*Industrial Quality Control & Material Microscopy*

*Metal Paint Electronics Coatings Ceramic Polymers Semiconductor Materials Gemology Museum Forensics*





## Use Video Port To Add 3D Imaging Capabilities

No objective or stage movement  
Sample space kept totally free  
No sample perturbation & vibrations

No motorization : maintenance free  
No microscope adaptation required  
Do not alter standard camera use



**ZeeScan on upright microscope  
with a Jenoptik Progres camera**




**ZeeScan on stereo microscope  
with a Micrometrics camera**



**ZeeScan on upright microscope  
with a fluorescence camera**

# MATERIAL MICROSCOPY


## High Resolution Digital Imaging



*Any c-mount cameras can be used with ZeeScan, with no alteration for the regular use of the camera. The ZeeScan optical assembly provides sharp & crisp digital for digital image documentation in high resolution.*

## Multiple Imaging Capabilities


**AutoFocus • Depth Measurement • Z-Stacking • Extended Depth of Field • 3D Reconstruction**



*ZeeScan performs all 3D critical tasks in material microscopy while using a standard upright or inverted microscope.*

## 3D Surface Metrology

**Surface Shape • Roughness • Waviness • Step Height**



*ZeeScan is the quickest and easiest way for precise surface topography measurements versus complex, bulky and expensive systems.*





## Software

- GetPhase® GUI software (included) is compatible with Windows 8, 7, XP & Vista. GetPhase® provides comprehensive tools from automatic acquisition to 2D / 3D image analysis, documentation and reports. An optional API / SDK includes ZeeScan acquisition controls, routines for Z-stack, 3D reconstruction, EDF, DIC, Phase, and 3D surface analysis.

### • Acquisition & Processing

- 2D / 3D Acquisition Wizard
- Auto Focus & Exposure
- Region-of-Interest
- Navigator
- Stitching
- Macro Recording

### • 2D/3D Display & Analysis

- BF, DF, Ph, DIC, 3D views
- Text & Graphics overlay
- 2D / 3D measurements
- Image fusion (EDF)
- Roughness ISO standards
- Step Height Measurements


### • Image Data Export & Report

- Project Archiving
- 3D Data in Excel Format
- 3D Data for 3<sup>rd</sup> Party Software
- Report Editor
- HTML Compatible Presentation

## Powerful Imaging Tool

Z-stacking of high resolution images can be automatically achieved providing image fusion image (Extended depth of Field image, Z depth measurement or 3D reconstruction. In addition, GetPhase provides 2D measurements and image documentation tools.


- ***Reveals finest structure details without specialized optics***
- ***On click Image documentation with multiple views***
- ***Automatic image fusion (Extended Depth of Field)***
- ***2D measurements & report***



## Fast & Accurate 3D Surface Metrology

ZeeScan with GetPhase performs 3D acquisition and analysis in a remarkable fast and easy way. Non contact optical surface profiling is highly repeatable.

- ***3D surface analysis in micrometer and nanometer range***
- ***Measurement capabilities from smooth to rough surfaces***
- ***ISO Roughness and step heights measurements***
- ***High throughput thanks to fast acquisition & processing time***





## Smart Hardware Architecture

No internal or external motorization, no additional accessories for the microscope are required, ZeeScan is connected to your PC using a single USB2 connection. Accurate calibration is achieved using an automated procedure and stored in an internal memory to prevent any losses.

Camera compatibility (not included)	Format 2/3" or Less, C-mount (see compatibility list)
Microscope Interface	Video Port – Recommended 1X C-mount adapter
PC Interface	USB 2.0
Power Supply	110 / 220 AC
Physical Dimensions (mm)	ZeeScan Head: 110(H) 80(W) 56(D) Control Unit: 40(H) 158(W) 150(D)
Weight	ZeeScan Head: 470g Control Unit: 150 g

## 3D Measurement Performance

Z range and resolution are objective and c-mount coupler magnification dependant. The table here under gives typical performance for standard objective magnification with 1X coupler. For any other magnification and /or c-mount coupler configuration, the following formulas can be applied:  $Z \text{ Range} = 16\text{mm} / (G_{\text{Obj}} * G_{\text{adapt}})^2$   $Z \text{ Resolution} = \text{Objective Depth Of Field} / 4$   
 $G_{\text{Obj}} = \text{Objective magnification}$   $G_{\text{adapt}} = \text{c-mount coupler magnification}$

Objective Mag / NA	Z Range (µm)	Z Resolution (µm)
5X / 0.10	640	0.32
10X / 0.25	160	0.08
20X / 0.45	40	0.02
50X / 0.8	6.4	0.003

*Z accuracy:* 1%

*Z Repeatability:* 0.35%

*Max slope:* 90°

*XY Spatial resolution determined by camera resolution and objective magnification*

## Roughness Measurement

12 analysis parameters are provided in total, including the frequently-used Ra (Sa), Rq (Sq), Rz (Sz), parameters. Parameters conform to ISO 4287, 25178 DIN 4768

Measuring range: Ra, Rq: 0.01-500µm Measuring accuracy:  $\leq \pm 10\%$  Repeatability:  $\leq 6\%$



2 Impasse de la Noisette - Bat B3, 91370 Verrières Le Buisson, France  
+33 9 54 03 05 43 • [contact@phaseview.com](mailto:contact@phaseview.com) • [www.phaseview.com](http://www.phaseview.com)